


More information available  
[www.hektros.com](http://www.hektros.com)



Please reach out to us for any questions about our products.



**Hennie Heidotting**  
[info@hektros.com](mailto:info@hektros.com)



Partners:



# TIPCUBE

## FULLY AUTOMATED PIPETTE TIP FILLING MACHINE

The TipCube is a compact pipette tip filling machine. It automatically racks 8 boxes in 20 min., through its innovative patented technology.

The TipCube is compatible with the most commonly used pipette tips of all brands and volumes, ranging from 0.1µl to 1,000µl.

Laboratories can save money, time and reduce plastic waste by using bulk pipette tips and reusable TipBoxes.



Save money



Save time




Less waste

## THE BEST PLACE FOR YOUR TIPS



The supplied TipBoxes from Hektros have a capacity of 96 Tips. They are available in red, yellow, and personalized colors.



If our boxes do not meet your requirements, we can supply you with the appropriate adapters for your 96 grid plate.

## SUPPORTED PIPETTE TIPS

The entire list is available on [www.hektros.com](http://www.hektros.com)



We prioritize the continual expansion and rigorous testing of our pipette tip portfolio to meet evolving customer needs. Here are our most popular selections:

BRAND	VOLUME	MODEL
Mettler Toledo	10 µl	Rainin RC UNV
SSI Bio	10 µl	SSI4137NSF
Eppendorf	200 µl	epT.I.P.S® Standard (yellow)
Gilson	200 µl	PIPETMAN DIAMOND Tips D200 EASYPACK
Inaoptica	200 µl	703C Microtips
Kirgen	200 µl	KG1232 Universal Fit Pipette Tips
SSI Bio	200 µl	SSI4230N00
Mettler Toledo	250 µl	Rainin LTS
Eppendorf	1,000 µl	epT.I.P.S® Standard

If you are not sure whether your pipette tips are supported, please feel free to contact us.

**<2** Minutes/Box  
FILLING TIME

**<2** Minutes  
SET UP TIME

**60** Watt  
MAX. POWER



## LET THE MACHINE DO A MACHINE'S JOB

The TipCube takes care of the tedious task of filling pipette tips, freeing up your employees to focus on their core responsibilities. With the TipCube, your team can concentrate on their actual work while the machine efficiently handles the filling process. This not only enhances productivity but also ensures that valuable time and effort are directed towards meaningful tasks, instead of mundanely filling boxes.

The TipCube offers an easy and intuitive operation. Its user-friendly design features a software interface that clearly explains each step, making it simple for anyone to operate. The rotating drum is mounted on a drawer and can be taken out, to fill it with bulk pipette tips. In the lower drawer, up to 8 TipBoxes or special Adapters can be inserted and filled within 20 minutes. Its compact design ensures it can fit seamlessly into different lab environments while maintaining high performance and reliability.

**0.1-1,000** µl  
VOLUME RANGE

**>1,000** Pipette tips  
DRUM CAPACITY